

A JAVID - CAROTID ARTERY CLAMP




A Carotid Artery Clamp
M077-100 Small
M077-105 Large



BLALOCK CLAMPS

B Blalock Clamp
M028-010 Small, 215mm
M028-015 Medium, 220mm
M028-020 Large, 225mm



BLALOCK - SERRATED JAW, CROSS ACTION

C Blalock Clamp, straight
G052-021 11mm jaw, 4cm
G052-023 14mm jaw, 4.5cm
G052-025 17mm jaw, 6cm
G052-027 27mm jaw, 7cm
G052-029 38mm jaw, 7.5mm
Blalock Clamp, curved
G052-020 11mm jaw, 4cm
G052-022 14mm jaw, 4.5cm
G052-024 17mm jaw, 6cm
G052-026 27mm jaw, 7cm
G052-028 38mm jaw, 7.5mm



SETTING FORCEPS - BULLDOG

D Setting Forceps - Bulldog
M064-675 23cm length



DENIS BROWNE CLAMP

E Denis Browne Clamp
M028-178 45° angle jaw
M028-180 90° angle jaw



HARKEN - BLOOD VESSEL CLAMP

F Harken - Blood Vessel Clamp 23cm
A - M028-301
B - M028-302
C - M028-303
D - M028-304



PRICE THOMAS BRONCHUS CLAMP, LONGITUDINAL SERRATIONS, WITH SPIKES

G Price Thomas Bronchus Clamp
G017-040 double angled, left 22cm
G017-042 double angled, right 22cm
G017-044 double angled, 22cm without spikes




SEMB - LIGATURE CLAMP

H Semb - Ligature Clamp
Distal eye for suture
M028-810 full curve, length 24cm
M028-812 shallow curve, length 24cm



SATINSKY - DOUBLE ANGLE



I Satinsky Clamps
Double angle, 26.5cm
M028-675 A, small angle
M028-680 B, large angle
Satinsky - Baby
Double angle
M028-685 length 15cm



A Debakey Clamps
M030-190 straight, 7.5cm
M030-192 straight, 8.5cm
M030-194 straight, 10.5cm
M030-196 straight, 12.0cm

M030-198 curved, 7.5cm
M030-200 curved, 8.5cm
M030-202 curved, 10.5cm
M030-204 curved, 12.0cm



Glover Clamps
M030-206 straight, 6cm
M030-208 straight, 7cm
M030-210 straight, 9cm
M030-212 straight, 11cm



M030-214 curved, 6cm
M030-216 curved, 7cm
M030-218 curved, 9cm
M030-220 curved, 11cm

B Cooley - Clamps
M030-018
 Anastomosis, 17cm

C Cooley - Derra Clamps
M030-020 small, 16cm
M030-022 medium, 17cm
M030-024 large, 18cm



D Cooley Multi Purpose Clamps
M030-040 30°, 17cm
M030-042 60°, 16cm
M030-044 90°, 15cm

E Cooley Anastomosis Clamps
M030-030
 double angle, 17cm
M030-032
 double angle, 17cm



F Cooley Ligature Clamps
M030-060 16cm
M030-070 17cm





G Cooley - Ductus Clamps
M030-050 straight, 17cm
M030-052 curved, 17cm
M030-054 angled, 17cm
M030-056 spoon shape, 17cm



H Cooley - Paediatric Clamp
M030-072 straight, 14cm
M030-074 angled, 14cm
M030-076 double angled, 14cm
M030-078 spoon shape, 16cm



DARDIK CLAMP
Multi Purpose





A M030-102 straight, 15cm

B M030-104 straight, 17cm



C M030-106 15° angle, 15cm

D M030-108 15° angle, 17cm





E M030-110 30° angle, 15cm

F M030-112 30° angle, 17cm









G M030-114 50° angle, 15cm

H M030-116 50° angle, 17cm



I M030-118 90° angle, 15cm

J M030-120 90° angle, 17cm



A Debakey Aorta Aneurysm Clamp
M030-150 full curve, 24.5cm
M030-152 medium curve
M030-154 slight curve, 26.5cm

B Debakey Aorta Aneurysm Clamp
M030-158 24cm
M030-160 30cm

C Debakey Aortic Clamp S - shape
M030-164 18.5cm
M030-166 30cm


D Debakey Aortic Clamp
M030-176 angled spoon, 26.5cm

E Debakey Aortic Clamp Shallow curve
M030-168 A, 24cm length
M030-170 B, 27cm length
M030-172 C, 31cm length
M030-174 D, 33.5cm length

F Debakey Derra Anastomosis Clamp
M030-182 A, 17cm
M030-184 B, 17cm
M030-186 C, 17cm

G Debakey Bainbridge Clamp
M030-690 straight, 40mm jaw, 15cm
M030-694 straight, 60mm jaw, 18cm
M030-692 curved, 40mm jaw, 15cm
M030-696 curved, 60mm jaw, 18cm


A DEBAKEY MULTI PURPOSE CLAMP



- A Debakey Multi Purpose Clamp**
Right Angle 90°
M030-221 A 18cm length
M030-222 B 19cm length
M030-223 C 23cm length
M030-224 D 27.5cm length

- Debakey Obtuse Angle 60°**
M030-226 A 20cm length
M030-227 B 22cm length
M030-228 C 24cm length
M030-229 D 30cm length

B DEBAKEY RING HANDLE BULLDOG CLAMP




- B Debakey - Ring Handle Bulldog Clamp**
M030-230 straight, 12cm
M030-232 curve, 11cm
M030-234 S-shape, 12cm
M030-236 angled, 13cm
M030-238 right angle, 10cm

- C Debakey - Curved Ligature Clamp**
M030-246 20cm length
M030-248 24cm length

- D Debakey - Shallow Curve Ligature Clamp**
M030-250 20cm length
M030-252 19cm length


- E Debakey - Full Curve Ligature Clamp**
M030-260 24cm length
M030-262 26.5cm length full serration

C DEBAKEY LIGATURE CLAMP



Curved

D DEBAKEY LIGATURE CLAMP



Shallow Curve


E DEBAKEY LIGATURE CLAMP



Full Curve

- F Debakey Coarctation Clamp**
M030-700 straight, 21cm
M030-702 straight, 23cm
M030-704 straight, 27cm
M030-706 straight, 30cm
M030-708 curved, 21cm
M030-710 curved, 23cm
M030-712 curved, 27cm
M030-714 curved, 30cm

F DEBAKEY COARCTATION CLAMP




G DEBAKEY SEMB LIGATURE CLAMP



- G Debakey Semb Ligature Clamp**
M030-264 25cm length

A DEBAKEY TANGENTIAL OCCLUSION CLAMP



A Debakey Tangential Occlusion Clamp

- M030-270 A** 20cm length
- M030-272 B** 25cm length
- M030-274 C** 26cm length
- M030-276 D** 27cm length
- M030-278 E** 28cm length

B GLOVER COARCTATION CLAMP



B Glover Coarctation Clamp

- M030-300** straight, 55mm jaw, 22cm
- M030-304** straight, 70mm jaw, 25cm
- M030-302** angled, 55mm jaw, 22cm
- M030-306** angled, 70mm jaw, 25cm
- M030-308** right angle curve, 18.5cm
- M030-312** shallow curve, 21cm
- M030-316** spoon shape, 21cm
- M030-310** right angle curve, 20cm
- M030-314** shallow curve, 24cm
- M030-318** spoon shape, 24cm


C Leland Jones Peripheral Clamp

- M030-675** straight, 19.5cm
- M030-680** angled, 15°, 19.5cm
- M030-685** angled, 30°, 19.5cm
- M030-687** angled, 45°, 18.5cm

C LELAND JONES PERIPHERAL CLAMP



D WEBER AORTIC FULL CURVE CLAMP



D Weber Aortic Full - curve Clamp

- M030-605** left, 26cm
- M030-610** right, 26cm


E Potts Coarctation Clamp

- M030-320** straight, half serration, 21cm
- M030-322** straight, full serration, 22cm
- M030-324** angled shank, half serration, 21cm
- M030-326** angled shank, full serration, 22cm

E POTTS COARCTATION CLAMP



F POTTS COARCTATION CLAMP




F Potts Coarctation Clamp

- M030-328** angled, half serration, 21cm
- M030-330** angled, full serration, 22cm

A SATINSKY CLAMP




B COOLEY NEONATAL - MINIATURE CLAMP




A Satinsky Clamp
 25.5cm length
M030-410 small
M030-415 medium
M030-420 large

B Cooley Neonatal Miniature Clamp
M030-501 No 1
M030-502 No 2
M030-503 No 3
M030-504 No 4
M030-505 No 5
M030-506 No 6
M030-507 No 7
M030-508 No 8
M030-509 No 9

C CASTANEDA NEONATAL - MULTI PURPOSE CLAMP




D MOSQUITO - ATRAUMATIC HEMOSTATIC FORCEPS



C Castaneda Neonatal Multi Purpose Clamp
M030-901 angled 30°, 13cm
M030-902 angled 50°, 13cm
M030-903 angled, 90°, 13cm
Anastomosis
M030-904 curved, small, 15cm
M030-905 curved, medium, 15cm
M030-906 curved, large, 15cm
M030-907 curved, extra large, 15cm

E PEAN - ATRAUMATIC HEMOSTATIC FORCEPS




D Mosquito - Atraumatic Hemostatic Forceps
M040-210 straight 12cm
M040-215 curved 12cm

E Pean - Atraumatic Hemostatic Forceps
M040-220 straight, 14cm
M040-222 straight, 16cm
M040-224 straight, 18cm
M040-226 straight, 20cm
M040-228 curved, 14cm
M040-230 curved, 16cm
M040-232 curved, 18cm
M040-234 curved, 20cm

F Debakey - Atraumatic Dissecting Forceps
G062-180 15cm length 1.5mm
G062-182 18cm length 1.5mm
G062-184 20cm length 1.5mm
G062-186 25cm length 1.5mm
G062-187 30cm length 1.5mm
G062-196 15cm length 2mm
G062-198 18cm length 2mm
G062-200 20cm length 2mm
G062-202 25cm length 2mm
G062-204 30cm length 2mm

G Debakey - Atraumatic Dissecting Forceps Angled
G062-204 15cm length 2mm
G062-206 18cm length 2mm
G062-208 20cm length 2mm
G062-210 25cm length 2mm
G062-212 30cm length 2mm

F DEBAKEY - ATRAUMATIC DISSECTING FORCEPS



G DEBAKEY - ATRAUMATIC DISSECTING FORCEPS ANGLED



1.5mm wide jaw, straight

2mm wide jaw, straight