







TUBE - LARYNGECTOMY



TUBE - LARYNGOSCOPY (EMERGENCY)



A Colledge Laryngectomy Tube
Silver, two tubes with shield & pilot

- E198-105** extra small
- E198-110** small
- E198-115** medium
- E198-120** large
- E198-125** extra large

B Colledge Laryngectomy Cage

E198-130 for use with Colledge Tubes

C Moure Lombards Laryngectomy Tube
Tube with one shield

- E198-205** 34fg
- E198-210** 38fg
- E198-215** 40fg

D Moure Lombards Laryngectomy Tube
Tube with two shields

- E198-250** 34fg
- E198-255** 36fg

E Shaw Colledge Composite Laryngectomy Tube

- | | |
|---|--|
| Tube with inner tube & pilot | Tube with two inner tubes & one pilot |
| E198-305 small | E198-350 small |
| E198-310 medium | E198-355 large |
| E198-315 large | |

F Cawthorne Laryngoscopy
Silver cannula, with blade, 10.2cm

E200-105

**REPAIR SERVICE
SURGICAL
INSTRUMENTATION**

Contact our
**CUSTOMER SERVICE
TEAM**

Tel : +44(0)1527 559010
Fax : +44(0)1527 559080
email : sales@judd-medical.co.uk

TRACHEOSTOMY TUBE REPAIR SERVICE

Comprehensive repair service for all makes of Tracheostomy Tubes.

Includes :


- Replacement of missing parts
- Total refurbishment
- Tubes returned in virtually new condition
- Rapid turnaround

Contact our
Customer Service Team
for further details


Tel: +44(0)1527 559010

Fax: +44(0)1527 559080

Email: sales@judd-medical.co.uk




TUBE - TRACHEOSTOMY

**A**
**Alder Hey
Tracheostomy Tubes**


Consisting of: Fenestrated outer tube with movable and lockable shield, fenestrated valved inner tube, plain inner tube, pilot and blocker. Silver

E202-012 12fg **E202-014** 14fg **E202-016** 16fg **E202-018** 18fg
E202-020 20fg **E202-022** 22fg **E202-024** 24fg

**B**
**Cubley
Tracheostomy Tubes**


Silver, tube with pilot

E202-316 16fg **E202-318** 18fg **E202-320** 20fg **E202-322** 22fg
E202-324 24fg **E202-326** 26fg **E202-328** 28fg **E202-330** 30fg
E202-332 32fg

**C**
**Edinburgh
Tracheostomy Tubes**


Silver, tube with pilot & moveable shield

E202-618 18fg **E202-620** 20fg **E202-622** 22fg **E202-624** 24fg
E202-626 26fg **E202-628** 28fg **E202-630** 30fg **E202-632** 32fg

**D**
**De Santi
Tracheostomy Tubes**


Silver, speaking aid valve

E202-424 24fg **E202-426** 26fg **E202-428** 28fg **E202-430** 30fg
E202-432 32fg **E202-434** 34fg **E202-436** 36fg **E202-438** 38fg

**E**
**Fuller
Tracheostomy Tubes**


Silver

E202-718 18fg **E202-720** 20fg **E202-722** 22fg **E202-724** 24fg

**F**
**Negus
Tracheostomy Tubes**


Silver, tube with pilot & one recessed & one plain valve

E202-814 14fg **E202-816** 16fg **E202-818** 18fg **E202-820** 20fg
E202-822 22fg **E202-824** 24fg **E202-826** 26fg **E202-828** 28fg
E202-830 30fg **E202-832** 32fg **E202-834** 34fg **E202-836** 36fg
E202-838 38fg **E202-840** 40fg

**G**
**Parker
Tracheostomy Tubes**

Silver, tube with pilot

E202-912 12fg **E202-914** 14fg **E202-916** 16fg **E202-918** 18fg
E202-920 20fg **E202-922** 22fg **E202-924** 24fg **E202-926** 26fg
E202-928 28fg **E202-930** 30fg **E202-932** 32fg **E202-934** 34fg

**H**
**Chevalier Jackson
Tracheostomy Tubes**

Tube with inner & pilot

Silver

E202-112 12fg **E202-114** 14fg **E202-116** 16fg **E202-118** 18fg
E202-120 20fg **E202-122** 22fg **E202-124** 24fg **E202-126** 26fg
E202-128 28fg **E202-130** 30fg **E202-132** 32fg **E202-134** 34fg
E202-136 36fg **E202-138** 38fg **E202-140** 40fg

TRACHEOSTOMY TUBE COLLARS

MARPAC



Correct fit and the patients comfort are assured with **Marpac** Tracheostomy Tube Collars. We offer a complete line to fit every size of patient - from “expand as you grow”, neonatal, paediatric to adult sizes.

Our easy to use collars are **made from soft foam** with **smooth sealed edges** to keep your patients comfortable. Plus they hold securely reducing the risk of accidental dislodgement or displacement of the tube and secondary complications.



TRACHEOSTOMY TUBE COLLARS

JM/TT/201 Adult, fits up to 21.5” neck **Box of 10**
JM/TT/203 Neonatal to Paediatric
JM/TT/204 Adult, fits up to 19” neck

JM/TT/201B Adult, fits up to 21.5” neck **Box of 100**
JM/TT/203B Neonatal to Paediatric
JM/TT/204B Adult, fits up to 19” neck

ENDOTRACHEAL TUBE TAPE



NEW ET Tape from **Marpac** takes less time to apply than conventional cloth tape, and **secures better**, which reduces the chances of the ET tube becoming dislodged or accidental extubation. The **skin friendly** 304 ET securement device includes zinc oxide, is **waterproof** and is **hypo-allergenic** for the patient’s comfort. ET Tape is individually packaged and has a unique **one piece design for quick application**. It’s made with a soft cotton/foam neck band and medical-grade adhesive tape. One end is enclosed in a rigid, disposable sleeve, which allows it to slide easily under the patient’s neck, and the scored backing is easy to peel off, even with protective gloves.

ENDOTRACHEAL TUBE TAPE

JM/ET/301 – Latex-free silk
JM/ET/304 – Latex-free cloth

Supplied in Boxes of 10

PAEDIATRIC ENDOTRACHEAL TUBE TAPE

JM/ET/305 – Latex-free silk
JM/ET/305B – Latex-free cloth
Size: 28½” long and 1” wide
 Supplied in Boxes of 10

