


## Cautery Burners

### PLATINUM TIPPED



**Point A**  
Code **F021-002** Straight, length 12.5cm  
Code **F021-102** Angled, length 12.5cm




**Point B**  
Code **F021-004** Straight, length 12.5cm  
Code **F021-104** Angled, length 12.5cm




**Point C**  
Code **F021-006** Straight, length 12.5cm  
Code **F021-106** Angled, length 12.5cm



**Point D**  
Code **F021-008** Straight, length 12.5cm  
Code **F021-108** Angled, length 12.5cm



**Point F**  
Code **F021-010** Straight, length 12.5cm  
Code **F021-110** Angled, length 12.5cm



**Point G**  
Code **F021-012** Straight, length 12.5cm  
Code **F021-112** Angled, length 12.5cm




**Point I**  
Code **F021-014** Straight, length 12.5cm  
Code **F021-114** Angled, length 12.5cm




**Point M**  
Code **F021-016** Straight, length 12.5cm  
Code **F021-116** Angled, length 12.5cm




**Point N**  
Code **F021-018** Straight, length 12.5cm  
Code **F021-118** Angled, length 12.5cm




**Point O**  
Code **F021-020** Straight, length 12.5cm  
Code **F021-120** Angled, length 12.5cm



**Straight**, length 12.5cm



**Angled**, length 12.5cm



**Point R**  
Code **F021-022** Straight, length 12.5cm  
Code **F021-122** Angled, length 12.5cm



**Point S**  
Code **F021-024** Straight, length 12.5cm  
Code **F021-124** Angled, length 12.5cm



**Point SP**  
Code **F021-026** Straight, length 12.5cm  
Code **F021-126** Angled, length 12.5cm



**Point LP**  
Code **F021-028** Straight, length 12.5cm  
Code **F021-128** Angled, length 12.5cm



**Point CP**  
Code **F021-030** Straight, length 12.5cm  
Code **F021-130** Angled, length 12.5cm

## Cautery Handles & Transformer


**Cautery Transformer** for mains power designed to modern electronic standards and conforming to EN 60601-1: EN 60601-1-2 1993 and are CE Marked.

- Finely graduated control
- Terminals to accept push-in connectors on cautery handle cable
- Convenient carrying handle
- 3m mains cable included (but not cable to handle)
- Size: 190mm x 220mm x 105mm high (including handle)
- Nett weight (including mains cable): 3.75kg

Code **JMR/650 Cautery Transformer**, 230V model, light duty for use with platinum tipped burners



### CAUTERY HANDLES




Mark Hovell

Code **E011-100**  
with press button switch

Pistol Grip

Code **E011-120** without illumination light & cords  
Code **E011-300** cords only, for cautery handles



Pistol Grip

Code **E011-115** fitted with illumination light, complete with integral lamp flex & cautery cord