

Gynae Theatre & Outpatients

January Product News

Ball Electrodes

Available **Sterile** & Non-sterile

Single Use
Fingerswitches
PMS250



Fingerswitch with 3 pin plug, 3m cable & blade electrode, **STERILE**
Boxed 25



PMS 1408
Pk of 3
110mm + ball 3mm
Sterile Code PMS 3408



PMS 1407
Pk of 3
110mm + ball 4mm
Sterile Code PMS 3407



PMS 1420
Pk of 3
110mm + ball 5mm
Sterile Code PMS 3407



PMS 1412
Pk of 3
110mm + ball 8mm
Sterile Code PMS 3407



PMS 1415
Pk of 3
150mm + ball 3mm
Sterile Code PMS 3415



PMS 1416
Pk of 3
150mm + ball 4mm
Sterile Code PMS 3416




PMS 1417
Pk of 3
150mm + ball 5mm
Sterile Code PMS 3417



PMS 1418
Pk of 3
150mm + ball 8mm
Sterile Code PMS 3418

Re-Usable
Colposcopy
Loops



See inside for
further details

Diathermy
Handcontrol
PMS 750

**SPECIAL
OFFERS**

Buy 10 Receive
20% DISCOUNT

Autoclavable
100 uses



View all our **SPECIAL OFFERS**
on our website @
www.judd-medical.co.uk/Special Offers.html

Re-Usable Colposcopy Loops



Safety guard: Any break in the outer black insulation is highlighted by the yellow insulation layer below.



Safety shield: A slip-on sleeve designed to prevent potential arcing between the electrode and fingerswitch.




Indicator grooves: Enables sterilisation records to be maintained using the appropriate systems.



Unique number: Enables sterilisation record to be made of the product. (up to 25 uses)


Packs of 5

1	PMS1509 8mm x 8mm, 11cm
2	PMS1510 10mm x 5mm, 11cm
3	PMS1511 10mm x 10mm, 11cm
4	PMS1512 12mm x 8mm, 11cm
5	PMS1501 12mm x 12mm, 11cm
6	PMS1513 15mm x 8mm, 11cm
7	PMS1502 15mm x 15mm, 11cm
8	PMS1514 20mm x 10mm, 11cm
9	PMS1503 20mm x 15mm, 11cm
10	PMS1515 25mm x 10mm, 11cm
11	PMS1517 25mm x 15mm, 11cm
12	PMS1504 25mm x 25mm, 11cm



Flexiline Pro

Single use Colposcopy Loop Electrodes for reduced adhesion and cleaner biopsies



The Prima Flexiline Pro Lletz system has been designed for increased precision in the treatment of cervical dysplasia. This is achieved by improving accuracy of removal of biopsy or lesion from the transformation zone.

Reinforced coated metal frame permits shaping of the active tip allowing the surgeon to manipulate the loop to the specific requirements of the patient, minimising removal of good tissue.

The symmetrical shape of the loop allows for clear vision of the active tip when in procedure, resulting in a more accurate tissue dissection, especially when working with larger loops where over removal can cause problems with misshaping of the cervix.

Non-stick ball electrodes minimise adhesion and prevent bleeders from resuming when concluding procedure. This further allows for fulguration to be achieved with greater success rates than with standard stainless steel ball electrodes.

Loop Electrodes


Packs of 10

PMS2610 10mm x 5mm, 11cm	PMS2614 20mm x 10mm, 11cm
PMS2600 8mm x 8mm, 11cm	PMS2602 15mm x 15mm, 11cm
PMS2611 10mm x 10mm, 11cm	PMS2615 25mm x 10mm, 11cm
PMS2612 12mm x 8mm, 11cm	PMS2603 20mm x 15mm, 11cm
PMS2613 15mm x 8mm, 11cm	PMS2617 25mm x 15mm, 11cm
PMS2601 12mm x 12mm, 11cm	PMS2604 25mm x 20mm, 11cm

Non-stick Ball Electrodes


Packs of 10

PMS2628 3mm ball electrode, 11cm
PMS2627 4mm ball electrode, 11cm
PMS2626 5mm ball electrode, 11cm



Colposcopy Loop Electrodes

Single Use



Packs of 25

1	PMS1740 Pink 10mm x 5mm loop
2	PMS1750 Aqua 10mm x 8mm loop
3	PMS1742 Green 10mm x 10mm loop
4	PMS1746 Orange 15mm x 8mm loop
5	PMS1743 Blue 15mm x 15mm loop
6	PMS1751 Grey 20mm x 12mm loop
7	PMS1752 Yellow 20mm x 13mm loop
8	PMS1745 Red 20mm x 15mm loop
9	PMS1744 Beige 20mm x 20mm loop
10	PMS1735 Light Green 5mm x 5mm square loop
11	PMS1731 Green 10mm x 4mm square loop
12	PMS1732 Brown 10mm x 8mm square loop

CUDA FIBEROPTICS

Judd Medical are **OFFICIAL UK DISTRIBUTORS** of the CUDA® range of products.

The CUDA® range consists of a quality, cost effective line of Fibre Optic Light Sources, Fibre Optic Head Lights and Cables. Our light sources are available with both Xenon and Metal Halide light options and provide solutions for Operating Theatres, Outpatient Departments and Clinics.

Our Headlights and Cables are fully compatible with all leading manufacturers light sources i.e. Olympus®, Storz®, Wolf® & BS/ACMI, through a unique system of interchangeable adaptors.

A full range of accessories is also available including; Replacement Elements/Bulbs, Floor stands, Head Light parts and Head Bands.

XENON LIGHT SOURCES

HEADLIGHTS

SPOTLIGHTS & ACCESSORIES

- XLS-300 - For use in Theatres and Suitable for Video Camera application
- XLS-180 - For use in Theatres and Out Patient Departments
- MLS-75T - For use in Out Patient Departments

- Deluxe Headlight with Ratchet
- Standard Headlight
- Deluxe Headlight with Clip

- 80mm & 110mm Straight
- 80mm & 110mm Bifurcated
- Floorstands
- Replacement bulbs/lamp elements



FIBRE LIGHT ADAPTORS

FIBRE LIGHT CABLES STRAIGHT OR BIFURCATED

- British Standard
- Wolf®/Dyonics®
- Karl Storz®
- Olympus®
- Storz®/Olympus®/GU®
- BS/ACMI



- British Standard
- Wolf®/Dyonics®
- Karl Storz®
- Olympus®
- Storz®/Olympus®/GU®
- BS/ACMI



Full CUDA® Fibreoptics catalogue available on request, or alternatively visit our website @ www.judd-medical.co.uk for further details.



Re-Usable & Single Use Electrosurgery Products

Judd Medical are delighted to announce that we have been appointed **EXCLUSIVE UK DISTRIBUTORS** for the Prima Medical range of Electrosurgery products. This partnership will allow Judd Medical to offer the same quality product, with the same price and product codes to the UK Market, with the advantage of large stocks held at our warehouse to ensure we continue to satisfy both existing Judd Medical and Prima Medical customers.

MONOPOLAR & BIPOLAR FORCEPS

COTSWOLD NEEDLES

UNIVERSAL BIPOLAR CABLES



EXCESS STOCK AVAILABLE!
Visit our website for more details



NEW PRODUCT


TWIN FOOTSWITCHES

DIATHERMY HANDCONTROLS

DISPERSIVE CABLES



SPECIAL OFFER
Visit our website for more details



SPECIAL OFFER
Visit our website for more details



Further details of the Prima range available on request, or alternatively visit our website @ www.judd-medical.co.uk for further details.

Judd Medical Limited,

Highfield House, 53 Worcester Road, Bromsgrove, Worcestershire, B61 7DN, UK

Tel: +44 (0)1527 559010 Fax: +44 (0)1527 559080

Email: sales@judd-medical.co.uk

www.judd-medical.co.uk