

**DEMONSTRATIONS  
AVAILABLE !**  
Judd Medical Limited  
are **EXCLUSIVE UK DISTRIBUTORS**  
for TAVA Surgical®  
Power Tools & Accessories



Also **NEW**  
from TAVA SURGICAL®  
**TOPS System**

*See inside for  
further details*



**Diathermy  
Handcontrol**  
PMS 750

**SPECIAL  
OFFERS**

Buy 10 Receive  
**20% DISCOUNT**

Autoclavable  
100 uses

Cable terminated with  
3-pin plug



View all our **SPECIAL OFFERS**  
on our website @  
[www.judd-medical.co.uk/Special Offers.html](http://www.judd-medical.co.uk/Special Offers.html)

**Power Tools & Accessories**  
from TAVA SURGICAL®

## New ..... Pulse Lavage

*The Tava Surgical® L31 Pulse Lavage System is an ergonomically designed and balanced, easy-to-use, pneumatic instrument that pumps fluids in distinct pulses at optimum pressure to remove blood and debris.*

*The L31 Pulse Lavage features an easy-to-load pump cartridge to facilitate easy operation.*

*A wide range of accessories including disposable tubing and nozzle sets designed for use during general purpose, hip and knee replacement surgeries are available.*

*Unique groove in the Pulse Lavage handle keeps the delivery tube neatly tucked in and out of the way.*

*Available with many popular hose connectors.*

### SPECIFICATIONS:

**Pulse Rate:** Variable to 1,700 pulses per minute by trigger and pressure adjustment

**Flow Rate:** Variable to .9 litres per minute

**Operating Pressure:** 75 to 110 psi (5-7bar)

*A full TAVA Surgical® catalogue is available on request, alternatively visit our website for full product details.*

Shown with the Knee Arthroplasty Set



+44(0)1527 559010



+44(0)1527 559080



sales@judd-medical.co.uk



www.judd-medical.co.uk

# New .....TOPS

## Tava Battery Operated Orthopaedic Power System

Every feature of Tava Surgical's Orthopaedic Power System is built with the surgeon in mind; from the ergonomically-designed handpiece to the quick-connect attachments, to the innovative sterile transfer system. Special patent pending features ensure that this system gives the surgeon smooth, precise and powerful operation and superior control.

This fully-functional, cordless device is equally suited to drilling, intramedullary reaming, sawing or inserting pins. It is tough enough and stable enough to go from cortical bone to cancellous bone and stay right on track. The Tava Surgical® Orthopaedic Power System offers an ergonomically designed handpiece, an all-inclusive, time-saving sterile transfer Power Module, and a wide variety of quick-connect, multipurpose attachments.

Judd Medical are able to supplement these drills with a full range of Burs, Blades, Pins, Wires & Twist Drills.

### Modular Handpiece



#### PERFORMANCE SPECIFICATIONS:

- Speed**
  - 100 to 700 rpm
- Torque Drill**
  - 40 in/lb (4.5Nm)
- Reamer**
  - 100 in/lb (11.3Nm)

### Sagittal/Oscillating Saw



#### PERFORMANCE SPECIFICATIONS:

- Speed**
  - 100 to 12,000 cpm
- Blade arc**
  - 5° arc

### Reciprocator/Sternal Saw



#### PERFORMANCE SPECIFICATIONS:

- Speed**
  - 100 to 17,000 cpm
- Blade Stroke**
  - 0.125"

### Modular Handpiece PM-X08-500

- Operates in forward and reverse with an added SAFE/OFF mode
- Reaches maximum speed in forward and reverse
- Couples to low profile, easy to use, reliable, quick-connect drilling and reaming attachments
- Fully cannulated to accept up to 0.156" / 5/32" / 4.0mm pins

### Sagittal/Oscillating Saw PM-X12-500

- Blade collet easy to use, quick install and release system
- 12 position rotating head at 30° intervals
- Innovative design - allows enhanced visibility for surgeon

### Reciprocator/Sternal Saw PM-X14-500

- Reciprocating Saw - Four 90° interval blade positions
- Sternum Saw - Two 180° interval blade positions

*The Tava Surgical® Orthopaedic Power System is specifically designed to address two of the most critical elements in orthopaedic surgery: reliable power and reduction of surgery time*



Other Products available from the Tava Surgical® Power Tools & Accessories range:

- **XK-PRO 100™ High Speed Drill System**
- **PneuMicro™ Wiredriver**
- **PneuMicro™ Small Bone Power System**

Full details on the TAVA Surgical® range is available on our website or alternatively contact our **CUSTOMER SERVICE TEAM**



SHUKLA™ Medical has over 160 years experience of providing world class products and services to the medical industry. Experienced engineers provide a full design and production facility ensuring their products are of the highest possible standard.

SHUKLA™ Medical have now partnered with Judd Medical Limited to provide **EXCLUSIVE UK DISTRIBUTION** of the SHUKLA™ range throughout the United Kingdom and together will continue to provide innovative products to assist clinical staff throughout all surgical specialities.

The SHUKLA™ Medical Xtract-All® products have been designed to offer a dedicated range of Implantation and Extraction Systems offering the Clinician a wide variety of solutions for specific surgical procedures.

**DEMONSTRATIONS AVAILABLE !**  
Contact our CUSTOMER SERVICE TEAM for further details +44(0)1527 559010

A full Shukla® catalogue is available on request, alternatively visit our website for full product details.

**Winquist III Universal IM Nail Extraction System**

Extraction of Intermedullary (IM) Nails including:-

- Cannulated Nails, Threaded and Non-Threaded
- Solid Nails
- Broken Nails (up to mid shaft)



**Universal Mini Ratcheting Screwdriver Set**


Removes & Installs small screws typically encountered in, Hand, Foot & Maxillofacial Surgery



Attachments include:  
■ 24 tips

**Universal Maxi Ratcheting Screwdriver Set**


Installs & extracts large bone screws



Attachments include:  
■ 8 Solid Drivers  
■ 8 Flexible Drivers

**Universal Broken, Stripped & Seized Screw Extraction System**

Removal of screws with stripped and damaged heads including screws cold fused on plates



Attachments include:  
■ 48 Extractors

# CUDA FIBEROPTICS

Judd Medical are **OFFICIAL UK DISTRIBUTORS** of the CUDA® range of products.

The CUDA® range consists of a quality, cost effective line of Fibre Optic Light Sources, Fibre Optic Head Lights and Cables. Our light sources are available with both Xenon and Metal Halide light options and provide solutions for Operating Theatres, Outpatient Departments and Clinics.

Our Headlights and Cables are fully compatible with all leading manufacturers light sources i.e. Olympus®, Storz®, Wolf® & BS/ACMI, through a unique system of interchangeable adaptors.

A full range of accessories is also available including; Replacement Elements/Bulbs, Floor stands, Head Light parts and Head Bands.

## XENON LIGHT SOURCES

## HEADLIGHTS

## SPOTLIGHTS & ACCESSORIES

- XLS-300 - For use in Theatres and Suitable for Video Camera application
- XLS-180 - For use in Theatres and Out Patient Departments
- MLS-75T - For use in Out Patient Departments

- Deluxe Headlight with Ratchet
- Standard Headlight
- Deluxe Headlight with Clip

- 80mm & 110mm Straight
- 80mm & 110mm Bifurcated
- Floorstands
- Replacement bulbs/lamp elements



## FIBRE LIGHT ADAPTORS

## FIBRE LIGHT CABLES STRAIGHT OR BIFURCATED

- British Standard
- Wolf®/Dyonics®
- Karl Storz®
- Olympus®
- Storz®/Olympus®/GU®
- BS/ACMI



- British Standard
- Wolf®/Dyonics®
- Karl Storz®
- Olympus®
- Storz®/Olympus®/GU®
- BS/ACMI



Full CUDA® Fibreoptics catalogue available on request, or alternatively visit our website @ [www.judd-medical.co.uk](http://www.judd-medical.co.uk) for further details.



## Re-Usable & Single Use Electrosurgery Products

Judd Medical are delighted to announce that we have been appointed **EXCLUSIVE UK DISTRIBUTORS** for the Prima Medical range of Electrosurgery products. This partnership will allow Judd Medical to offer the same quality product, with the same price and product codes to the UK Market, with the advantage of large stocks held at our warehouse to ensure we continue to satisfy both existing Judd Medical and Prima Medical customers.

### MONOPOLAR & BIPOLAR FORCEPS

### UNIVERSAL MONOPOLAR CABLES

### UNIVERSAL BIPOLAR CABLES




**EXCESS STOCK AVAILABLE!**  
Visit our website for more details

**NEW PRODUCTS**

### TWIN FOOTSWITCHES

### DIATHERMY HANDCONTROLS

### DISPERSIVE CABLES



**SPECIAL OFFER**  
Visit our website for more details

**SPECIAL OFFER**  
Visit our website for more details

Further details of the Prima range available on request, or alternatively visit our website @ [www.judd-medical.co.uk](http://www.judd-medical.co.uk) for further details.