


Gynae Theatre & OPD

March Product News

Flexiline Pro

Single use Colposcopy Loop Electrodes for reduced adhesion and cleaner biopsies



The Prima Flexiline Pro Lletz system has been designed for increased precision in the treatment of cervical dysplasia. This is achieved by improving accuracy of removal of biopsy or lesion from the transformation zone.

Reinforced coated metal frame permits shaping of the active tip allowing the surgeon to manipulate the loop to the specific requirements of the patient, minimising removal of good tissue.

The symmetrical shape of the loop allows for clear vision of the active tip when in procedure, resulting in a more accurate tissue dissection, especially when working with larger loops where over removal can cause problems with misshaping of the cervix.

Non-stick ball electrodes minimise adhesion and prevent bleeders from resuming when concluding procedure. This further allows for fulguration to be achieved with greater success rates than with standard stainless steel electrodes.

Single use Non-Stick Ball Electrodes



Non-stick Ball Electrodes	
Packs of 10	
PMS2628	3mm ball electrode, 11cm
PMS2627	4mm ball electrode, 11cm
PMS2626	5mm ball electrode, 11cm

Loop Electrodes	
Packs of 10	
PMS2610	10mm x 5mm, 11cm
PMS2600	8mm x 8mm, 11cm
PMS2611	10mm x 10mm, 11cm
PMS2612	12mm x 8mm, 11cm
PMS2613	15mm x 8mm, 11cm
PMS2601	12mm x 12mm, 11cm

PMS2614	20mm x 10mm, 11cm
PMS2602	15mm x 15mm, 11cm
PMS2615	25mm x 10mm, 11cm
PMS2603	20mm x 15mm, 11cm
PMS2617	25mm x 15mm, 11cm
PMS2604	25mm x 20mm, 11cm

PMS2614

PMS2602

PMS2615

All electrodes 11cm long

NEW UNIVERSAL Monopolar Cable

One Cable Fits all and a choice of colours



Safety sleeve - see the breaks before they occur

- Double insulated cable
- Built-in cable restraints
- Highly flexible cable/connector

Contact our CUSTOMER SERVICE TEAM for full details

TWIN FOOTSWITCHES **SPECIAL OFFER**




Part No: FS355FQ

CALL US FOR SPECIAL OFFER PRICES

View all our **SPECIAL OFFERS** on our website @ www.judd-medical.co.uk/Special Offers.html

Colposcopy Loop Electrodes Single Use




Packs of 25

1	PMS1740 Pink 10mm x 5mm loop
2	PMS1750 Aqua 10mm x 8mm loop
3	PMS1742 Green 10mm x 10mm loop
4	PMS1746 Orange 15mm x 8mm loop
5	PMS1743 Blue 15mm x 15mm loop
6	PMS1751 Grey 20mm x 12mm loop
7	PMS1752 Yellow 20mm x 13mm loop
8	PMS1745 Red 20mm x 15mm loop
9	PMS1744 Beige 20mm x 20mm loop
10	PMS1735 Light Green 5mm x 5mm square loop
11	PMS1731 Green 10mm x 4mm square loop
12	PMS1732 Brown 10mm x 8mm square loop

Ball Electrodes Single Use - STERILE

1	PMS3410 3mm Ball, 5cm	Boxed 25
2	PMS3403 4mm Ball, 5cm	Boxed 25
3	PMS3435 5mm Ball, 5cm	Boxed 25
4	PMS3405 8mm Ball, 5cm	Boxed 25
5	PMS3408 3mm Ball, 11cm	Boxed 25
6	PMS3407 4mm Ball, 11cm	Boxed 25
7	PMS3420 5mm Ball, 11cm	Boxed 25
8	PMS3412 8mm Ball, 11cm	Boxed 25
9	PMS3415 3mm Ball, 15cm	Boxed 10
10	PMS3416 4mm Ball, 15cm	Boxed 10
11	PMS3417 5mm Ball, 15cm	Boxed 10
12	PMS3418 8mm Ball, 15cm	Boxed 10




Diathermy Quivers


Diathermy Quiver
4cm diameter x 23cm length
G022-905



Diathermy Quiver
6cm diameter x 23cm length
G022-910



Julian Taylor Diathermy Quiver
with drainage holes,
4.0cm diameter x 23cm length
G022-915




Diathermy Quiver
with removable screw bottom
G022-970
5cm diameter x 20cm length
G022-975
6cm diameter x 23cm length
G022-985
6cm diameter x 37.5cm length



Diathermy Quiver
6cm diameter
x 37.5cm length
G022-950



Diathermy Quiver
Telescopic for laparoscopic
instruments,
6cm diameter x 40cm length
G022-920




TOWEL CLIPS

Ball & Socket Towel Clip
G090-290 11cm



Other types & sizes available

Fingerswitch with 3 pin
plug, 3m cable & blade
electrode, **STERILE**
Boxed 25



Single Use
Fingerswitches
PMS250





Judd Medical Limited are OFFICIAL UK DISTRIBUTORS for CUDA® Fibreoptics

SPOTLIGHTS & ACCESSORIES

- 80mm & 110mm Straight
- 80mm & 110mm Bifurcated
- Floorstands
- Replacement bulbs/lamp elements



HEADLIGHTS

- Deluxe Headlight with Ratchet
- Standard Headlight
- Deluxe Headlight with Clip elements



XENON LIGHT SOURCES

- XLS-300Q - For use in Theatres and Suitable for Video Camera application
- XLS-180Q - For use in Theatres and Out Patient Departments
- MLS-75T - For use in Out Patient Departments



FIBRE LIGHT ADAPTORS

- British Standard
- Wolf®/Dyonics®
- Karl Storz®
- Olympus®
- Storz®/Olympus®/GU®
- BS/ACMI



FIBRE LIGHT CABLES STRAIGHT OR BIFURCATED

- British Standard
- Wolf®/Dyonics®
- Karl Storz®
- Olympus®
- Storz®/Olympus®/GU®
- BS/ACMI



CUDA® FIBRE LIGHT CABLES CALL US FOR SPECIAL OFFER PRICES

Full CUDA® Fibreoptics catalogue available on request, or alternatively visit our website @ www.judd-medical.co.uk for further details.



Re-Usable & Single Use Electrosurgery Products

Judd Medical are delighted to announce that we have been appointed **EXCLUSIVE UK DISTRIBUTORS** for the Prima Medical range of Electrosurgery products. This partnership will allow Judd Medical to offer the same quality product, with the same price and product codes to the UK Market, with the advantage of large stocks held at our warehouse to ensure we continue to satisfy both existing Judd Medical and Prima Medical customers.

MONOPOLAR & BIPOLAR FORCEPS



EXCESS STOCK AVAILABLE!
Visit our website for more details

COTSWOLD NEEDLES



UNIVERSAL BIPOLAR CABLES



NEW PRODUCT

DIATHERMY HANDCONTROLS


PMS 750Q Buy 10 Receive 20% DISCOUNT

Autoclavable 100 uses



Cable terminated with 3-pin plug

COLPOSCOPY LOOPS



DISPERSIVE CABLES



Further details of the Prima range available on request, or alternatively visit our website @ www.judd-medical.co.uk for further details.

Judd Medical Limited,
Highfield House, 53 Worcester Road, Bromsgrove, Worcestershire, B61 7DN, UK
Tel: +44 (0)1527 559010 Fax: +44 (0)1527 559080
Email: sales@judd-medical.co.uk

www.judd-medical.co.uk