

### Replacement Burs & Blades

Judd Medical offer an extensive range of **REPLACEMENT BURS & BLADES**, featured are a small selection of our most popular. For further details on the range, please contact our CUSTOMER SERVICE TEAM Tel: +44(0)1527 559010

CE Marked  
Blades supplied Sterile for Single Use  
**DISCOUNTS** are available for Contract/Call off Orders

CAT NUMBER    HEAD DIAMETER    FLUTES    LENGTH

#### ZIMMER® - Type Long Burs



**JZB-330**    5.0 mm    8    95.0mm  
Carbide

#### ZIMMER® - Type Extra Long Burs




**JZB-236**    4.0mm    8    70.0mm  
Carbide

A full copy of our Replacement Bur and Blade catalogue is available on request


CAT NUMBER    BLADE DEPTH    BLADE WIDTH    BLADE THICKNESS

#### STRYKER® - Type System 6 Sagittal Saw Blades



**ST2108-100**    73.0 mm    25.0 mm    0.64 mm

#### STRYKER® - Type Oscillating Blades




**JSO-652**    81.5 mm    12.7mm    1.20 mm

#### HALL/ZIMMER® - Type Large Reciprocating Blades



**JZR-279**    68.6 mm    .60mm    12.5 mm

#### HALL/ZIMMER® MICRO-CHOICE® - Type Flat Sagittal Blades



**JZS-237**    25.5 mm    14.0mm    0.40 mm

**FREE Headlight**  
with every XLS-300Q  
& XLS-180Q  
Light Source



**TWIN FOOTSWITCHES**    **SPECIAL OFFER**



**CALL US FOR SPECIAL OFFER PRICES**

View all our **SPECIAL OFFERS** on our website @ [www.judd-medical.co.uk/Special Offers.html](http://www.judd-medical.co.uk/Special Offers.html)

**Tava Surgical® is a division of Brasseler USA™ Inc.**  
**Using research engineers and consulting physicians, Tava Surgical® have developed a line of Surgical Power Tools & Accessories for the Operating Theatre.**

**Judd Medical Limited have now partnered with Tava Surgical® to provide EXCLUSIVE DISTRIBUTION in the UK of the TAVA Surgical® range of Power Tools & Accessories.**

## Pulse Lavage System & Accessories

The L31 Pulse Lavage is an easy-to-use, pneumatic instrument that pumps fluids in distinct pulses at optimum pressure to remove blood and debris.



## DEMONSTRATIONS AVAILABLE !

Contact our CUSTOMER SERVICE TEAM for further details +44(0)1527 559010

## XK-PRO 100™ High Speed Drill System

The XK-PRO 100™ High Speed Drill System is a compact, lightweight design that offers new technological and performance standards for a combination of surgical procedures. Rotating at 100,000 rpm's, the pneumatic powered high-speed drill system exemplifies qualities to include smooth, precise cutting, superior fingertip control, minimal operational noise level, and is virtually oil free.



## PneuMicro™ Small Bone Power System

The Small Bone Power System offers precision with ease of operation. The system can be utilized for all small bone procedures to include orthopaedic, oral maxillofacial, podiatric and otology techniques.

RECIPROCATING SAW



OSCILLATING SAW



SAGITTAL SAW



## TOPS - Battery Operated Orthopaedic Power System

Every feature of the TOPS is designed with the surgeon in mind; from the ergonomically designed hand-piece to the quick-connect attachments, to the innovative sterile transfer system. This fully-functional, cordless device is equally suited to drilling, intramedullary reaming, sawing or inserting pins.



## PneuMicro™ Wiredriver

The PneuMicro™ Wiredriver from Tava Surgical® brings dependable lube-free technology to the surgeon. This pneumatically-powered unit, with a comfort grip handle design, places precise, variable speed control at your fingertips.





SHUKLA™ Medical has over 160 years experience of providing world class products and services to the medical industry. Experienced engineers provide a full design and production facility ensuring their products are of the highest possible standard.

SHUKLA™ Medical have now partnered with Judd Medical Limited to provide **EXCLUSIVE UK DISTRIBUTION** of the SHUKLA™ range throughout the United Kingdom and together will continue to provide innovative products to assist clinical staff throughout all surgical specialities.

The SHUKLA™ Medical Xtract-All® products have been designed to offer a dedicated range of Implantation and Extraction Systems offering the Clinician a wide variety of solutions for specific surgical procedures.

## DEMONSTRATIONS AVAILABLE !

Contact our CUSTOMER SERVICE TEAM for further details +44(0)1527 559010

A full SHUKLA™ catalogue is available on request, alternatively visit our website for full product details.

### Winquist III Universal IM Nail Extraction System

Extraction of Intermedullary (IM) Nails including:-

- Cannulated Nails, Threaded and Non-Threaded
- Solid Nails
- Broken Nails (up to mid shaft)



### Universal Mini Ratcheting Screwdriver Set

Removes & Installs small screws typically encountered in, Hand, Foot & Maxillofacial Surgery



Attachments include:  
■ 24 tips

### Universal Maxi Ratcheting Screwdriver Set


Installs & extracts large bone screws



Attachments include:  
■ 8 Solid Drivers  
■ 8 Flexible Drivers

### Universal Broken, Stripped & Seized Screw Extraction System

Removal of screws with stripped and damaged heads including screws cold fused on plates



Attachments include:  
■ 48 Extractors

# CUDA<sup>®</sup> FIBEROPTICS

Judd Medical are **OFFICIAL UK DISTRIBUTORS** of the CUDA<sup>®</sup> range of products.

The CUDA<sup>®</sup> range consists of a quality, cost effective line of Fibre Optic Light Sources, Fibre Optic Head Lights and Cables. Our light sources are available with both Xenon and Metal Halide light options and provide solutions for Operating Theatres, Outpatient Departments and Clinics.

Our Headlights and Cables are fully compatible with all leading manufacturers light sources i.e. Olympus<sup>®</sup>, Storz<sup>®</sup>, Wolf<sup>®</sup> & BS/ACMI, through a unique system of interchangeable adaptors.

A full range of accessories is also available including; Replacement Elements/Bulbs, Floor stands, Head Light parts and Head Bands.

## XENON LIGHT SOURCES

## HEADLIGHTS

## SPOTLIGHTS & ACCESSORIES

- XLS-300Q - For use in Theatres and Suitable for Video Camera application
- XLS-180Q - For use in Theatres and Out Patient Departments
- MLS-75T - For use in Out Patient Departments

- Deluxe Headlight with Ratchet
- Standard Headlight
- Deluxe Headlight with Clip

- 80mm & 110mm Straight
- 80mm & 110mm Bifurcated
- Floorstands
- Replacement bulbs/lamp elements



**CUDA<sup>®</sup> FIBRE LIGHT CABLES**  
CONTACT US FOR SPECIAL OFFER PRICES

## FIBRE LIGHT ADAPTORS

## FIBRE LIGHT CABLES STRAIGHT OR BIFURCATED

- British Standard
- Wolf<sup>®</sup>/Dyonics<sup>®</sup>
- Karl Storz<sup>®</sup>
- Olympus<sup>®</sup>
- Storz<sup>®</sup>/Olympus<sup>®</sup>/GU<sup>®</sup>
- BS/ACMI



- British Standard
- Wolf<sup>®</sup>/Dyonics<sup>®</sup>
- Karl Storz<sup>®</sup>
- Olympus<sup>®</sup>
- Storz<sup>®</sup>/Olympus<sup>®</sup>/GU<sup>®</sup>
- BS/ACMI

Full CUDA<sup>®</sup> Fibreoptics catalogue available on request, or alternatively visit our website @ [www.judd-medical.co.uk](http://www.judd-medical.co.uk) for further details.



## Re-Usable & Single Use Electrosurgery Products

Judd Medical are delighted to announce that we have been appointed **EXCLUSIVE UK DISTRIBUTORS** for the Prima Medical range of Electrosurgery products. This partnership will allow Judd Medical to offer the same quality product, with the same price and product codes to the UK Market, with the advantage of large stocks held at our warehouse to ensure we continue to satisfy both existing Judd Medical and Prima Medical customers.


### MONOPOLAR & BIPOLAR FORCEPS

### UNIVERSAL MONOPOLAR CABLES

### BLADE ELECTRODES



**EXCESS STOCK AVAILABLE!**  
Visit our website for more details



**NEW PRODUCTS**

### FINGERSWITCHES

### UNIVERSAL BIPOLAR CABLES

### DISPERSIVE CABLES

#### PMS250

Fingerswitch with 3 pin plug, 3m cable & blade electrode, **STERILE** Boxed 25



Further details of the Prima range available on request, or alternatively visit our website @ [www.judd-medical.co.uk](http://www.judd-medical.co.uk) for further details.

## Judd Medical Limited,

Highfield House, 53 Worcester Road, Bromsgrove, Worcestershire, B61 7DN, UK

Tel: +44 (0)1527 559010 Fax: +44 (0)1527 559080

Email: [sales@judd-medical.co.uk](mailto:sales@judd-medical.co.uk)

[www.judd-medical.co.uk](http://www.judd-medical.co.uk)