

TRAUMA SURGERY

SHUKLA™ Medical - Screwdriver Implantation & Extraction Systems MINI RATCHETING SCREWDRIVER SYSTEM

SHUKLA™ Medical has over 160 years experience of providing world class products and services to the medical industry. Experienced engineers provide a full design and production facility ensuring their products are of the highest possible standard.

SHUKLA™ Medical have now partnered with Judd Medical Limited to provide **Exclusive Distribution** of the SHUKLA™ range throughout the United Kingdom and together will continue to provide innovative products to assist Clinical Staff throughout all surgical specialities.

The SHUKLA™ Medical Xtract-all® products have been designed to offer a dedicated range of **Implantation and Extraction Systems** offering the Clinicians a wide variety of solutions for specific surgical procedures.

Installs & Removes:

- Small Screws typically encountered in:
 - Hand Surgery
 - Foot Surgery
 - Maxillofacial Surgery

Attachments include:

- 24 Tips
 - Osteotome
 - 3 Torx® Drives
 - Phillips® Drive
 - Clutch Drive
 - 4 Centre Drives
 - 2 Slotted Drives
 - 3 Cruciform Drives
 - 3 Hex Drives
 - 4 Stripped Screw Extractors
 - 2 Broken Screw Extractors
 - Ratcheting Handle



NEW PRODUCTS
Judd Medical are **EXCLUSIVE UK DISTRIBUTORS**

WINQUIST III UNIVERSAL IM NAIL EXTRACTION SYSTEM

UNIVERSAL BROKEN, STRIPPED & SEIZED SCREW EXTRACTION SYSTEM

Removes:

- Cannulated Nails
 - Threaded
 - Non Threaded
- Solid Nails
- Broken Nails (up to mid-shaft)




Attachments include:

- 24 Extractors
 - 51 Single Use Extractor Tips
 - 5 Multi Use Extractors
- Bone Chisel
- Screwdriver Set
 - Including patented Reverse Captive Twist™ Heads
 - Straight & T Ratcheting Handles
 - 18 Blades
- 5lb. Mallet & Striking Frame

DEMONSTRATIONS ARE AVAILABLE ON THE FULL RANGE OF SHUKLA™ PRODUCTS
Contact our **CUSTOMER SERVICE TEAM** to book your appointment

Removes:

- M2 to M8 Size Broken & Stripped Screws
- Cold Weld (fused) Stripped Screws
Capable of applying over 250 in-lbs of Torque



















Attachments include:

- 48 Extractors
- 3 Extensions
- Breaker Bar
- Jacobs Chuck Attachment
- Ratcheting T-Handle

CALL or FAX your order to our Customer Service Team !!!!

Tel: +44(0)1527 559010 Fax: +44(0)1527 559080


<p>MCINDOE DISSECTING FORCEP serrated, 15cm G062-340</p> 	<p>ARTHUR (STIEGLITZ) SPLINTER FORCEP straight, 14cm G082-102</p> 	<p>HALSTEAD MOSQUITO ARTERY FORCEP curved, 12.5cm G052-106</p> 	<p>RAMPLEY SPONGE HOLDER straight, 24cm G084-224</p> 
<p>IRIS SCISSORS sh/sh, straight, 11.5cm G140-270</p> 	<p>IRIS SCISSORS sh/sh, curved, 11.5cm G140-275</p> 	<p>KILNER NEEDLE HOLDER 13.5cm G124-112</p> 	<p>GILLIES SKIN HOOK round handle, 17cm G104-110 small G104-120 large</p> 
<p>KILNER CATSPAW (SENN-MILLER) standard, double ended G134-540</p> 	<p>CUSCO SPECULUM G200-020 small G200-025 medium G200-035 large G200-015 virgin</p> 	<p>BLOCK END DISSECTING FORCEP 12.5cm G062-040</p> 	<p>SIMS SPECULUM G200-070 small G200-075 medium G200-080 large</p> 
<p>WEST RETRACTOR self retaining, 3:4 sharp prongs, 14cm G136-720</p> 	<p>ALLIS TISSUE FORCEP 2:3 teeth, 15cm G088-023</p> 	<p>HARTMAN AURAL FORCEP crocodile action, standard, serrated jaw, satin finish, 7.0mm x 1.25mm E036-415</p> 	<p>HARTMAN AURAL FORCEP crocodile action, extra fine, serrated jaw, satin finish, 4.0mm x 0.8mm E036-420</p> 

A & E Product Mailing

<p>TILLEY DRESSING FORCEP serrated, angled, 16.5cm E054-105</p> 	<p>LISTER SINUS FORCEP serrated jaw, 12.5cm E078-105</p> 	<p>THUDICUM NASAL SPECULUM E182-501 No 1 12mm E182-502 No 2 14mm E182-503 No 3 16mm E182-504 No 4 18mm</p> 	<p>GRAEFE IRIS FORCEP 1:2 teeth, 1mm straight, 7cm H122-118</p> 
<p>TUMARKIN AURAL SPECULUM Available sizes : No 0 - No 7</p> 	<p>MINI FRIEDMAN BONE NIBBLING FORCEP slight curve, 12.5cm E068-103</p> 	<p>CARTRIDGE SYRINGE ward quality, 2.2cc P320-100M</p> 	<p>SPENCER WELLS CURVED ARTERY FORCEP G052-374 12.5cm G052-394 18cm</p> 
<p>PLIER WITH SIDE CUTTER cuts up to 1.6mm, soft wires, 15cm T800-045</p> 	<p>STERILE PIN COVERS 10-1-001 Yellow Up to 1.6mm 10-1-002 Blue Up to 4.0mm 10-1-003 Red Up to 6.0mm</p> 	<p>KALUSKAR FISHBONE FORCEP 12.5cm E061-750</p> 	<p>MAYO NEEDLE HOLDER 25cm G124-156</p> 
<p>LISTER BANDAGE SCISSORS 18cm S/S, CE Marked, autoclavable JMS/MC1 Pk of 5</p> 	<p>MULTI-CUT SCISSORS S/S blade, CE Marked, autoclavable, washable, 18cm JMS/MC1 Pk of 10</p> <p>Available:- Red, Green, Yellow & Blue</p> 	<p>RING CUTTER G018-105 16.5cm G018-110 Spare blade</p> 	<p>TROUSSEAU TRACHEA DILATING FORCEP E056-205 Child, 10.5cm E056-210 Adult, 14cm</p> 

For more information on the range of products supplied by **JUDD MEDICAL LTD** please go to **www.judd-medical.co.uk** or call us on **+44(0)1527 559010**

BALL & SOCKET TOWEL CLIP
G090-290 11cm
G090-295 14cm



TREVES DISSECTING FORCEP
 1:2 teeth, 12.5cm
G062-484



GILLIES DISSECTING FORCEP
 1:2 teeth, 15cm
G062-268



DISPOSABLE CURETTES
 Safe Ear Currettes
 All types
 Colour coded, Box 50



MAYO INSTRUMENT PINS
JM/IP10 - 10cm
JM/IP15 - 15cm




MAGILL INTRODUCING FORCEP
A090-050 Adult, 23cm
A090-100 Child, 19cm



DEBAKEY DISSECTING FORCEP
 1.5mm tip, 16cm, T/C jaw
G063-180




JOBSON HORNE PROBE
 double ended, 14cm
E150-305




BAG CLIPS
 stainless steel,
 25mm x 50mm
G090-100 Pk 10




SLING CLIPS
 stainless steel,
 65mm jaw
G090-150 Pk 10



MAGNETIC MATS
 Available Sterile & Reusable
80855 - 12" x 8"
80885 - 12" x 16"



MACHS
 To assist with stabilising
 leads & tubes
 Colour coded to suit your
 requirements



Issue Oct 09 JD



Judd Medical Limited,
 Highfield House, 53 Worcester Road, Bromsgrove, Worcestershire, B61 7DN, UK
 Tel: +44 (0)1527 559010 Fax: +44 (0)1527 559080
 Email: sales@judd-medical.co.uk

www.judd-medical.co.uk

