

TRACHEOSTOMY TUBE COLLARS

Judd Medical Limited
Official UK Distributors

- **Easy to fit**
- **Comfortable**
- **Secure**



TRACHEOSTOMY TUBE COLLARS

JM/TT/201	Adult, fits up to 21.5" neck	Box of 10
JM/TT/203	Neonatal to Paediatric	
JM/TT/204	Adult, fits up to 19" neck	
JM/TT/201B	Adult, fits up to 21.5" neck	Box of 100
JM/TT/203B	Neonatal to Paediatric	
JM/TT/204B	Adult, fits up to 19" neck	

Correct fit and the patients comfort are assured with Marpac Tracheostomy Tube Collars. We offer a complete line to fit every size of patient - from "expand as you grow", neonatal, paediatric to adult sizes.

Our easy to use collars are **made from soft foam** with **smooth sealed edges** to keep your patients comfortable. Plus they hold securely reducing the risk of accidental dislodgement or displacement of the tube and secondary complications.

Bariatric Tracheostomy Collars - 'Custom Made Sizes to Order'

Contact our CUSTOMER SERVICE TEAM for more information

Judd Medical Limited, Highfield House, 53 Worcester Road, Bromsgrove, Worcestershire, B61 7DN, UK

ENDOTRACHEAL TUBE TAPE

NEW PRODUCT

PAEDIATRIC ENDOTRACHEAL TUBE TAPE

Quick, effective securement for your Paediatric Patients



NEW ET Tape from **Marpac** takes less time to apply than conventional cloth tape, and **secures better**, which reduces the chances of the ET tube becoming dislodged or accidental extubation.

The **skin friendly** ET securement device includes zinc oxide, is **waterproof** and is **hypo-allergenic** for the patient's comfort. ET Tape is individually packaged and has a unique **one piece design for quick application**. It's made with a soft cotton/foam neck band and medical-grade adhesive tape. One end is enclosed in a rigid, disposable sleeve, which allows it to slide easily under the patient's neck, and the scored backing is easy to peel off, even with protective gloves.

PAEDIATRIC ENDOTRACHEAL TUBE TAPE

Size: 28½" long and 1" wide

JM/ET/305 – Latex-free
Supplied in Boxes of 10

JM/ET/305B
Supplied in Boxes of 10



ENDOTRACHEAL TUBE TAPE



- **Quick to fit**
- **Effective**
- **Comfortable**

ENDOTRACHEAL TUBE TAPE

JM/ET/301 - Latex-free silk
JM/ET/304 - Latex-free cloth

Supplied in Boxes of 10

