

New Range of Orthopaedic Sagittal Saw Blades

Judd Medical Limited is proud to introduce its **NEW RANGE** of precision Replacement Saw Blades suitable for use in Orthopaedic procedures.

Judd Medical has been a manufacturer and distributor of quality saw blades and burs since 1997 and have established itself as a leader in the UK Market by combining product quality and performance with cost effective pricing structures.

We have an extensive portfolio of Burs and Blades that are compatible with all leading manufacturers power tool systems. All our products are supplied in sterile packaging and are sterilised via Gamma Irradiation offering a five year sterile shelf life complete with batch traceability labelling.


For further information, or to request a full Bur and Blade Catalogue please contact our **Customer Service Team** on **+44(0)1527 559010**

JUDD REF	Usable Length mm	Blade Width mm	Blade Thickness mm	Cut Thickness mm	Tooth Pitch mm (TPI)	De Soutter Ref
----------	------------------	----------------	--------------------	------------------	----------------------	----------------


Illustrations not to scale

Sagittal Saw Blades


NEW RANGE




J10683	25	19	0.4	0.7	1.1 (23)	S83-106
--------	----	----	-----	-----	----------	---------



J10883	25	15	0.4	0.7	1.1 (23)	S83-108
--------	----	----	-----	-----	----------	---------




J10983	48	15	0.4	0.7	1.15 (22)	S83-109
--------	----	----	-----	-----	-----------	---------




J15183	25	10	0.7	1.0	1.1 (23)	S83-151
--------	----	----	-----	-----	----------	---------




J15383	25	6.4	0.7	1.0	1.1 (23)	S83-153
--------	----	-----	-----	-----	----------	---------



J20080	90	25	0.8	0.97	2.1 (12)	S80-200
--------	----	----	-----	------	----------	---------



J20680	90	25	1.20	1.27	2.1 (12)	S80-206
--------	----	----	------	------	----------	---------



J20780	90	25	1.30	1.47	2.1 (12)	S80-207
--------	----	----	------	------	----------	---------


ORTHOPAEDIC SAGITTAL SAW BLADES

**NEW
RANGE**


JUDD REF	Usable Length mm	Blade Width mm	Blade Thickness mm	Cut Thickness mm	Tooth Pitch mm (TPI)	De Soutter Ref
----------	------------------	----------------	--------------------	------------------	----------------------	----------------

Illustrations not to scale


Sagittal Saw Blades




J10280	70	25	0.8	1.0	1.7 (15)	S80-102
--------	----	----	-----	-----	----------	---------




J10080	87	25	0.8	1.0	1.7 (15)	S80-100
--------	----	----	-----	-----	----------	---------



J10380	60	12	0.6	0.8	1.7 (15)	S80-103
--------	----	----	-----	-----	----------	---------




J10480	87	19	0.8	1.0	1.7 (15)	S80-104
--------	----	----	-----	-----	----------	---------




J10980	32	12	0.5	0.8	1.3 (20)	S80-109
--------	----	----	-----	-----	----------	---------




J10082	45	12	0.5	0.9	1.6 (16)	S82-100
--------	----	----	-----	-----	----------	---------



J10682	95	25	0.81	1.0	1.6 (16)	S82-106
--------	----	----	------	-----	----------	---------



J00299	95	25	1.20	1.27	2.1 (12)	S99-002
--------	----	----	------	------	----------	---------



J00399	95	19	1.20	1.27	2.1 (12)	S99-003
--------	----	----	------	------	----------	---------


ORTHOPAEDIC SAGITTAL SAW BLADES

NEW RANGE


JUDD REF	Usable Length mm	Blade Width mm	Blade Thickness mm	Cut Thickness mm	Tooth Pitch mm (TPI)	De Soutter Ref
----------	------------------	----------------	--------------------	------------------	----------------------	----------------

Illustrations not to scale


Sagittal Saw Blades



J10680	87	25	1.20	1.27	2.2 (11)	S80-106
--------	----	----	------	------	----------	---------



J82099	90	25.4	1.20	1.27	1.7 (15)	S99-820
--------	----	------	------	------	----------	---------



J82599	80	25.4	1.20	1.27	1.7 (15)	S99-825
--------	----	------	------	------	----------	---------

Along with the new products, Judd Medical continue to offer its extensive range of Burs & Blades, a selection of which are detailed below.

JUDD REF	Usable Length mm	Blade Width mm	Blade Thickness mm	Cut Thickness mm	Tooth Pitch mm (TPI)	Equivalent Manu Ref	De Soutter Ref
----------	------------------	----------------	--------------------	------------------	----------------------	---------------------	----------------

Illustrations not to scale

Micro Reciprocating Saw Blade

Saw Blades




JZR-021	39	9.4	0.36	0.6	0.8 (32)	5053-021/221 Zimmer®	S22-021
---------	----	-----	------	-----	----------	----------------------	---------

Reciprocating Saw Blades



JZR-061	89	-	1.35	1.7	1.4 (18)	5052-061/261 Zimmer®	S22-061
---------	----	---	------	-----	----------	----------------------	---------



JSR-275	78	11	0.8	1.0	1.6 (16)	277-96-275 Stryker®	S57-275
---------	----	----	-----	-----	----------	---------------------	---------




JSR-325	78	11	0.8	1.0	1.6 (16)	277-96-325 Stryker®	S57-325
---------	----	----	-----	-----	----------	---------------------	---------

ORTHOPAEDIC SAW BLADES


Illustrations not to scale

Small Oscillating Saw Blade



JZS-032


Fine Tooth Micro Sagittal Saw Blade




JZS-038

JUDD REF	Usable Length mm	Blade Width mm	Blade Thickness mm	Cut Thickness mm	Tooth Pitch mm (TPI)	Equivalent Manu Ref	De Soutter Ref
JZS-032	25.4	6	0.4	0.65	1.3 (20)	5053-032/232 Zimmer®	S17-032
JZS-038	25	9.4	0.4	0.6	0.7 (36)	5053-038/238 Zimmer®	S17-038

Micro Sagittal Saw Blades



JSP-125



JSP-414

JSP-125	31	9	0.38	0.45	0.8 (32)	2296-3-125 Stryker®	S59-125
JSP-414	25	5.5	0.38	0.45	0.8 (32)	2296-3-414 Stryker®	S59-414

Sagittal Saw Blades



JMX-821	100	25	1.12	1.27	-	4125-127-100 Stryker®	S60-406
---------	-----	----	------	------	---	-----------------------	---------




JMX-823	90	25	1.12	1.27	-	4125-127-090 Stryker®	S60-403
---------	----	----	------	------	---	-----------------------	---------



JSO-619	95	25	1.25	1.27	2.1 (12)	2108-218 Stryker®	S60-218
---------	----	----	------	------	----------	-------------------	---------



JSO-682	85	19	1.25	1.27	2.1 (12)	2108-182 Stryker®	S60-182
---------	----	----	------	------	----------	-------------------	---------




JSO-683	90	29	1.25	1.27	2.1 (12)	2108-183 Stryker®	S60-183
---------	----	----	------	------	----------	-------------------	---------

Oscillating Saw Blades




JMO-124	50	30	0.6	1.1	2.0 (12)	L124 3M®	S32-127
---------	----	----	-----	-----	----------	----------	---------




JMO-131	70	29	1.0	1.1	1.3 (20)	L131 3M®	S32-131
---------	----	----	-----	-----	----------	----------	---------

Sagittal Saw Blades




JMS-134



JMS-136

JMS-134	32	10.4	0.4	0.7	1.0 (26)	K134 3M®	S33-134
JMS-136	20	7.5	0.4	0.7	0.9 (29)	K136 3M®	S33-136



J1200-005



J1400-001

J1200-005	19.8	8.9	0.3	0.5	0.7 (36)	1200-005 MicroAire®	S72-005
J1400-001	19.8	12.7	0.3	0.5	0.7 (36)	1200-001 MicroAire®	S72-001

Issue July 09 JD



Judd Medical Ltd,
Highfield House, 53 Worcester Road,
Bromsgrove, Worcestershire, B61 7DN, UK
Tel: +44 (0)1527 559010 Fax: +44 (0)1527 559080
Email: sales@judd-medical.co.uk

www.judd-medical.co.uk

