

The 'Piano Pedal' Footswitch


PLASTIC COATED ALUMINIUM ELECTROSURGERY FOOTSWITCH

- Non-slip operation
- Robust
- Ideal for use in the sitting position
- Easy to clean
- Non - corrosive plastic finish



'Piano Pedal' Footswitch

PLASTIC COATED ALUMINIUM ELECTROSURGERY FOOTSWITCH



Outer guard directs foot contact



Built in cable restraint



Integral carrying handle



Ideal for use with all theatre footwear

Ordering Details - all cables 5 meters

For use with the following manufacturer's electrosurgery generators




Valleylab	FS3206
Bard / Davol	FS3041
Conmed / Aspen	FS3250

Eschmann 300 & 400 Series	
	FS3207

Eschmann 800 Series	
	FS3307

Erbe	
	FS3720

No Connector	
	FS3355F