


*New, from Tava Surgical...*

# **PneuMicro™** Wiredriver

*For Precise, Dependable Control and Power*



The PneuMicro™ Wiredriver from Tava Surgical brings dependable lube-free technology to the surgeon. This pneumatically-powered unit, with a comfort grip handle design, places precise, variable speed control at your fingertips.

The PneuMicro™ Wiredriver's motor and valve delivers the control and power necessary for drilling and wire fixation, and our collet design offers expanded versatility of drilling with quick connect drills.

## *An Array of New, Outstanding Features:*

- Fatigue-reducing, ergonomic handle design.
- Lightweight and compact.
- Highly responsive, variable speed trigger.
- Improved, robust, high-performance motor design.
- Lube-free motor design for reduced maintenance.
- Expanded wire/pin diameter range...from 0.7 mm to 2.0 mm.
- Convenient, single hand directional control switch.
- Improved valve design for enhanced directional control with "positive touch."

**Judd**  
medical

## Accessories



Wire Guard



Wire Brush



Cleaning Cap

## Wiredriver Specifications

<b>Catalog Number:</b>	PM-M11-200
<b>Speed:</b>	2,100 rpm
<b>Torque:</b>	45 in-oz. (31.8 N-cm)
<b>Wire/Pin Diameter Range:</b>	0.028 in - 0.078 in (0.7 mm - 2.0 mm)
<b>Operating Pressure:</b>	90 - 110 psi (6.2 - 7.6 kg/cm <sup>2</sup> )

## Product Ordering Information

<b>Handpiece</b>	PM-M11-200*	Wiredriver
<i>* Purchase of handpiece includes Wire Guard and Wire Brush</i>		
<b>Accessories</b>	KM-A90-200	Air/Nitrogen Hose (Hall®) - 10 ft
	KM-A90-201	Cleaning Cap
	PM-M11-201	Wire Guard
	PM-M11-202	Wire Brush

## Other Tava PneuMicro™ Products

### Handpieces

PM-M09-200	High Speed Drill
PM-M10-200	Medium Speed Drill
PM-M12-200	Sagittal Saw
PM-M13-200	Oscillating Saw
PM-M14-200	Reciprocating Saw

### Attachments

PM-M10-901	Medium Bur Guard
PM-M10-902	Long Bur Guard
PM-M10-903	Extra Long Bur Guard

### Accessories

KM-A90-200-15	Air/Nitrogen Hose (Hall®) - 15 ft
KM-A91-000	Foot Control
KM-A95-000	Sterilization Case - Small (11 in x 7 in x 1 in)
MI-121-001	Adapter - DISS to Schrader
PM-A95-000	Sterilization Case - Large (12 in x 7 in x 4 in)
PM-M12-001	Blade Wrench

Judd Medical Limited  
Highfield House, 53 Worcester Road, Bromsgrove, Worcestershire,  
B61 7DN, UK  
Tel: +44(0)1527 559010 Fax: +44(0)1527 559080  
Email: sales@judd-medical.co.uk

**Judd**  
medical