

**Prima**  
Medical Limited

ELECTROSURGERY  
**COLPOSCOPY**



[www.judd-medical.co.uk](http://www.judd-medical.co.uk)



## Designed and Manufactured by Prima Medical Ltd


All rights reserved. No part of this publication may be reproduced or transmitted in any form or by any other means, whether electronic or mechanical including by photocopy, recording or information storage and retrieval system, without permission in writing from the publisher. Authorisation to photocopy items for internal or personal use, or the internal or personal use of specific clients is granted by Prima Medical Ltd.

### All Products are Latex Free

Copyright © 2006, Prima Medical Ltd.

## CONTENTS

Biopsy Punches & Gynae Hook	2
Loop Electrodes (Re-Usable)	3
Fingerswitches & Ball Electrodes (Re-Usable)	4
Loop Electrodes (Single Use)	5
Fingerswitches & Ball Electrodes (Single Use)	6



[www.judd-medical.co.uk](http://www.judd-medical.co.uk)



**BIOPSY PUNCH**


1 Biopsy Punch Tishler-Morgan Re-usable pack 1 PMS1000

**GYNAE HOOK**

2 Gynae Hook 250mm Single Use pack 50 PMS1790



## LOOPS




## RE-USABLE COLPOSCOPY ELECTRODES

1	Loop Electrode, 8 mm × 8 mm loop, 11cm,	pack 5	PMS1509
2	Loop Electrode, 10 mm × 5 mm loop, 11cm,	pack 5	PMS1510
3	Loop Electrode, 10 mm × 10 mm loop, 11cm,	pack 5	PMS1511
4	Loop Electrode, 12 mm × 8 mm loop, 11cm,	pack 5	PMS1512

5	Loop Electrode, 12 mm × 12 mm loop, 11cm,	pack 5	PMS1501
6	Loop Electrode, 15 mm × 8 mm loop, 11cm,	pack 5	PMS1513
7	Loop Electrode, 15 mm × 15 mm loop, 11cm,	pack 5	PMS1502
8	Loop Electrode, 20 mm × 10 mm loop, 11cm,	pack 5	PMS1514

9	Loop Electrode, 20 mm × 15 mm loop, 11cm,	pack 5	PMS1503
10	Loop Electrode, 25 mm × 10 mm loop, 11cm,	pack 5	PMS1515
11	Loop Electrode, 25 mm × 15 mm loop, 11cm,	pack 5	PMS1517
12	Loop Electrode, 25 mm × 25 mm loop, 11cm,	pack 5	PMS1504

Size may vary by 5%



## PRIMA SAFETY FEATURES

## Safety guard



Any break in the outer black insulation is highlighted by the yellow insulation layer below

## PRIMA SAFETY FEATURES

## Safety shield



A slip-on sleeve designed to prevent potential arcing between the electrode and fingerswitch.

## PRIMA SAFETY FEATURES

## Indicator grooves



Enables sterilisation record to be made on the product. (up to 25 uses)

## UNIQUE NUMBER




Enables sterilisation records to be maintained using the appropriate systems.


## FINGERSWITCH



## BALL ELECTRODES



## LOOP ELECTRODES




## RE-USABLE FINGERSWITCH

- 1 100 Plus Fingerswitch, 3m cable terminated with 3 pin plug,  
(Autoclavable for a maximum of 100 cycles) pack 5 PMS750

## RE-USABLE BALL ELECTRODES

- 2 3 mm Ball Electrode, 11 cm, pack 3 PMS1408  
3 5 mm Ball Electrode, 11 cm, pack 3 PMS1420





LOOPS



SINGLE USE COLPOSCOPY ELECTRODES


Size may vary by 5%

1 Loop Electrode, 10mm x 5mm loop, Pink, pack 25	PMS1740	5 Loop Electrode, 15mm x 15mm loop, Blue, pack 25	PMS1743	9 Loop Electrode, 20mm x 20mm loop, Beige, pack 25	PMS1744
2 Loop Electrode, 10mm x 8mm loop, Aqua, pack 25	PMS1750	6 Loop Electrode, 20mm x 12mm loop, Grey, pack 25	PMS1751	10 Loop Electrode, 5mm x 5mm square loop, Light Green, pack 25	PMS1735
3 Loop Electrode, 10 mm x 10mm loop, Green, pack 25	PMS1742	7 Loop Electrode, 20mm x 13mm loop, Yellow, pack 25	PMS1752	11 Loop Electrode, 10mm x 4mm square loop, Green, pack 25	PMS1731
4 Loop Electrode, 15mm x 8mm loop, Orange, pack 25	PMS1746	8 Loop Electrode, 20mm x 15mm loop, Red, pack 25	PMS1745	12 Loop Electrode, 10mm x 8mm square loop, Brown, pack 25	PMS1732





#### BALL ELECTRODES



#### SINGLE USE FINGERSWITCH

1 Fingerswitch with 3 pin plug, 3m cable, boxed 25 PMS255

#### SINGLE USE BALL ELECTRODES

2 Ball Electrode, 3 mm, 12cm pack 25 PMS1761

3 Ball Electrode, 5 mm, 12cm pack 25 PMS1762



[www.judd-medical.co.uk](http://www.judd-medical.co.uk)

Call today  
to receive your  
**FREE**  
Prima catalogue

