

## The non-blocking orthopaedic suction system

- Non-blocking with easy specimen removal
- Extra large filtration chamber
- Two reversible tips



# Freeflo

THE NON-BLOCKING SUCTION SYSTEM

- **Non blocking**
- **Easy specimen removal**
- **Two reversible tips 15cm / 23cm**

It has been acknowledged that when using suction systems during orthopaedic surgery, problems occur.

When drawing blood product away from an open area, traditional suction units may become blocked due to blood coagulation and polymerisation, typically between the handle and suction bottle.

The Freeflo handle draws fluid at high speed, trapping debris inside the filter chamber, ensuring that this does not block the system elsewhere with all the subsequent inconvenience.

## Extra large filtration chamber

The unique filter handle, designed for the collection of debris before it causes blocking, is fifty per-cent larger in capacity than previous designs.

## Ordering Details -



**200-102**

Freeflo non-blocking suction unit with 15cm and 23cm tips (box of 10)



**200-202**

Freeflo non-blocking suction unit with 15cm and 23cm tips and 3m suction tube (box of 10)



15cm tip

23cm tip