

Orthoclean

orthopaedic brushes for the preparation of bone surfaces prior to implantation



• soft grip handles

• nylon bristles –
to scour bone surface

• spirally wound bristles –
to lift debris

Code	Description	Pack
OBR100	Femoral brush for power tools/Jacobs chuck, sterile 27cm	10
OBR150	Femoral brush with handle, sterile, 34cm	10
OBR200	Acetabular brush for power tools/Jacobs chuck, sterile 19cm	10
OBR250	Acetabular brush with handle, sterile 24cm	10
OBR300	Femoral brush for power tools/Jacobs chuck, sterile 27cm	100
OBR350	Femoral brush with handle, sterile 34cm	100
OBR400	Acetabular brush for power tools/Jacobs chuck, sterile 19cm	100
OBR450	Acetabular brush with handle, sterile 24cm	100